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M-6FN DYE

【MODEL】

M-6FN

【PACKAGE SPECIFICATIONS】

48mL×1,12mL×4

【PRINCIPIO】

This product is formulated to dye the NRBCs and leukocytes with the assist of M-6LN LYSE. It participates in the measurement of NRBC-related and basophil-related parameters which is accomplished by using the laser scatter and the SF Cube cell analysis technology (3D analysis using information from scatter of laser light at two angles and fluorescence signals).

【INTENDED USE】

The M-6FN DYE participates in the measurement of Baso-related and NRBC-related parameters together with M-6LN LYSE.

【ACTIVE INGREDIENTS】

Fluorochrome..... ≤0.01%
Ethylene Glycol..... 94-99 %

【PRODUCT SPECIFICATIONS】

1. Appearance: Blue and transparent liquid without deposits, suspended grains or flocks.
2. Blank count results: The blank count results of the reagent tested on the Mindray Auto Hematology Analyzer shall meet the requirements of Table 1.

Table 1 Blank Count Requirements

Parameter	Blank Count Requirements
WBC-N	≤0.2×10 ⁹ /L

【STORAGE AND EXPIRATION DATE】

This product shall be stored in ventilated room free from corrosive gas and light at 2°C–30°C with the humidity lower than 90%. The shelf life is 1 year.
The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 60 days.

【APPLICABLE INSTRUMENT】

This product applies to
BC-6000/BC-6000Plus/BC-6100/BC-6100Plus/BC-6200/BC-6600Plus/BC-6700Plus/
BC-6800Plus/BC-7500[N]/BC-7500[NR]/BC-7500[N] CRP/BC-7500[NR] CRP/
BC-7300[N] CRP/BC-7300[NR] CRP/BC-7500[N] CS/BC-7500[NR] CS Auto Hematology Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

【INSTRUCTIONS FOR USE】

1)Restore the product to usage temperature;
2) Open the external packing and insert the pickup tube into the reagent pack based on the matching between the reagent pack label and the inlet label of the analyzer cap assembly connector;
3) Screw the cap assembly tight and replace the reagent according to the operator's manual of the analyzer;
4) Run a blank count and check the results. If the results meet the blank count requirements defined in the operator's manual of the analyzer, the newly installed reagent can be used for sample analysis See the operator's manual of the analyzer for details.

【WARNINGS】

Do not inhale and/or ingest. In case of inhalation or ingestion, induce vomiting and seek medical advice immediately.
Avoid skin contact. In case of skin contact, flush the affected area with plenty of water immediately.
Avoid eye contact. In case of eye contact, flush the affected area with plenty of water and seek medical advice immediately.

【NOTICE】

1) This product shall only be used on matched Mindray Auto Hematology Analyzer. Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
2) If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one. Dispose of the waste, residual and polluted packing based on local regulations.
3) After the product is opened, it shall be connected to the cap assemblies at once and shall not be left open in the air.
4) Keep the product away from light and high temperature.
5) The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
6) This product must be operated by clinical professionals, physician or laboratorian trained by Mindray or Mindray authorized personnel.

【INTELLECTUAL PROPERTY STATEMENT】

SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD. (hereinafter called Mindray) owns the intellectual property rights to this Mindray product and this manual. This manual may refer to information protected by copyright or patents and does not convey any license under the patent rights or copyright of Mindray, or of others.
Mindray intends to maintain the contents of this manual as confidential information. Disclosure of the information in this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.
Release, amendment, reproduction, distribution, rental, adaptation, translation or any other derivative work of this manual in any manner whatsoever without the written permission of Mindray is strictly forbidden.

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Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address: Mindray Building, Keji 12th Road South, High-tech industrial park, Nanshan, Shenzhen 518057, P.R.China
E-mail Address: service@mindray.com
Tel: +86 755 81888998
Fax: +86 755 26582680
EC-Representative: Shanghai International Holding Corp. GmbH(Europa)
Address: Eiffestraße 80, 20537 Hamburg, Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

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